



# **Keep Going**



Applied Home Healthcare Equipment, Westlake, OH 44145




### **Keeps Your Patients**

active!

Whether they are out on the town or home with family—patients can always have their oxygen with them. Restoring freedom and control of their day is a big part of what makes the OxyGo™ so special.

Patients can power an OxyGo™ by plugging it into a wall outlet, DC outlet in their car, or for hours at a time with each rechargeable battery. Weighing in at less than 5 pounds, and with an extraordinary battery life of 9 hours, the OxyGo™ is completely wearable.



## OxyGo™ Comparison Chart

	Lightweight							
	OxyGo™ single cell battery	OxyGo™ double cell battery	Respironics SimplyGo Mini® std. battery	Respironics SimplyGo Mini® ext. battery	Inova Labs LifeChoice Activox 4L®	Inova Labs Activox Pro®	CAIRE <sup>®</sup> FreeStyle <sup>™</sup>	CAIRE® Focus™
O <sub>2</sub> Capacity	1050 ml/min	1050 ml/min	1000 ml/min	1000 ml/min	450 ml/min	450 ml/min	498 ml/min	332 ml/min
Overnight	~	<b>&gt;</b>	×	×	<b>&gt;</b>	<b>~</b>	×	×
Single Solution	~	<b>~</b>	×	×	×	×	×	×
User Replaceable Batteries	~	~	~	~	×	×	×	*
Intelligent Delivery Technology: Flow Settings: 1,2,3,4,5	~	*	×	×	×	×	×	×
Weight	4.8 lbs.	5.8 lbs.	5.0 lbs.	6.0 lbs.	4.8 lbs.	4.3 lbs.	4.4 lbs.	2.3 lbs.
	2.2 kg.	2.6 kg.	2.3 kg.	2.7 kg.	2.0 kg.	1.95 kg.	2.0 kg.	1.0 kg.
Battery Duration	UP TO 4.5 hours	UP TO 9.5 hours	UP TO 4.5 hours	UP TO 9 hours	UP TO 8.25 hours	UP TO 7.25 hours	UP TO 3.5 hours	UP TO 3.0 hours
Sound Level at setting 2	39 dB.	39 dB.	43 dB.	43 dB.	44 dB.	46 dB.	44 dB.	Unspecified

View the clinical studies: bitly.com/OxyGoClinical

#### **OxyGo**<sup>™</sup> Portable Oxygen Concentrator

Give us a call to learn more about how OxyGo<sup>™</sup> can keep your patients going, and keep your costs down!

888-327-7301 or med@applied-inc.com

# 1400-1000

# 1400-1000-16 (system with double cell battery)



#### www.oxygo.life

